



MRGPRX2 Antibody, HRP conjugated

Product Code	CSB-PA839361LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB1
Immunogen	Recombinant Human Mas-related G-protein coupled receptor member X2 protein (286-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mas-related G-protein coupled receptor member X2, MRGPRX2, MRGX2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MRGPRX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.