



C2orf50 Antibody

Product Code	CSB-PA839365LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LR7
Immunogen	Recombinant Human Uncharacterized protein C2orf50 protein (1-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C2orf50, C2orf50
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C2orf50