




# CFAP52 Antibody

<b>Product Code</b>	CSB-PA839781LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8N1V2
<b>Immunogen</b>	Recombinant Human Cilia- and flagella-associated protein 52 protein (101-400AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cilia- and flagella-associated protein 52 (WD repeat-containing protein 16) (WD40-repeat protein up-regulated in HCC), CFAP52, WDR16 WDRPUH
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cancer
<b>Target Names</b>	CFAP52

## Image



### Western blot

All lanes: CFAP52 antibody at 4µg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: K562 whole cell lysate

Lane 3: U937 whole cell lysate


Lane 4: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 69, 70, 61 kDa

Observed band size: 69 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA839781LA01HU at dilution of 1:100

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.