



PRSS35 Antibody, Biotin conjugated

Product Code	CSB-PA839792LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N3Z0
Immunogen	Recombinant Human Inactive serine protease 35 protein (201-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive serine protease 35, PRSS35, C6orf158
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRSS35
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.