



METTL25 Antibody, HRP conjugated

Product Code	CSB-PA839813LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6Q8
Immunogen	Recombinant Human Methyltransferase-like protein 25 protein (206-455AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methyltransferase-like protein 25 (EC 2.1.1.-), METTL25, C12orf26
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	METTL25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.