





USHBP1 Antibody

Product Code	CSB-PA839815LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6Y0
Immunogen	Recombinant Human Usher syndrome type-1C protein-binding protein 1 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Usher syndrome type-1C protein-binding protein 1 (USH1C-binding protein 1) (AIE-75-binding protein) (MCC-2) (Mutated in colon cancer protein 2), USHBP1, AIEBP MCC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USHBP1


Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA839815LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA839815LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA839815LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)