



USHBP1 Antibody, Biotin conjugated

Product Code	CSB-PA839815LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6Y0
Immunogen	Recombinant Human Usher syndrome type-1C protein-binding protein 1 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Usher syndrome type-1C protein-binding protein 1 (USH1C-binding protein 1) (AIE-75-binding protein) (MCC-2) (Mutated in colon cancer protein 2), USHBP1, AIEBP MCC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USHBP1