



USHBP1 Antibody, Biotin conjugated

Product Code	CSB-PA839815LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6Y0
Immunogen	Recombinant Human Usher syndrome type-1C protein-binding protein 1 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Usher syndrome type-1C protein-binding protein 1 (USH1C-binding protein 1) (AIE-75-binding protein) (MCC-2) (Mutated in colon cancer protein 2), USHBP1, AIEBP MCC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USHBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.