



OR5L1 Antibody, FITC conjugated

Product Code	CSB-PA839894LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGL2
Immunogen	Recombinant Human Olfactory receptor 5L1 protein (161-196AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 5L1 (OST262) (Olfactory receptor OR11-151), OR5L1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	OR5L1