



USP33 Antibody, HRP conjugated

Product Code	CSB-PA840583LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEY7
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 33 protein (251-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 33 (EC 3.4.19.12) (Deubiquitinating enzyme 33) (Ubiquitin thioesterase 33) (Ubiquitin-specific-processing protease 33) (VHL-interacting deubiquitinating enzyme 1) (hVDU1), USP33, KIAA1097 VDU1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	USP33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.