



DEFB104A Antibody, HRP conjugated

Product Code	CSB-PA840974LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTQ1
Immunogen	Recombinant Human Beta-defensin 104 protein (23-72AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 104 (Beta-defensin 4) (BD-4) (DEFB-4) (hBD-4) (Defensin, beta 104), DEFB104A; DEFB104B, DEFB104 DEFB4;
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	DEFB104A