



KLK6 Antibody, HRP conjugated

Product Code	CSB-PA842165LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92876
Immunogen	Recombinant Human Kallikrein-6 protein (22-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-6 (EC 3.4.21.-) (Neurosin) (Protease M) (SP59) (Serine protease 18) (Serine protease 9) (Zyme), KLK6, PRSS18 PRSS9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KLK6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.