



KLK6 Antibody, Biotin conjugated

Product Code	CSB-PA842165LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92876
Immunogen	Recombinant Human Kallikrein-6 protein (22-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-6 (EC 3.4.21.-) (Neurosin) (Protease M) (SP59) (Serine protease 18) (Serine protease 9) (Zyme), KLK6, PRSS18 PRSS9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KLK6