



BORCS7 Antibody, Biotin conjugated

Product Code	CSB-PA842639LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96B45
Immunogen	Recombinant Human BLOC-1-related complex subunit 7 protein (2-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BLOC-1-related complex subunit 7 (Diaskedin), BORCS7, C10orf32
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BORCS7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.