



ARL8A Antibody

Product Code	CSB-PA842642LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BM9
Immunogen	Recombinant Human ADP-ribosylation factor-like protein 8A protein (20-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-like protein 8A (ADP-ribosylation factor-like protein 10B) (Novel small G protein indispensable for equal chromosome segregation 2), ARL8A, ARL10B GIE2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ARL8A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.