



ARL8A Antibody

Product Code	CSB-PA842642LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BM9
Immunogen	Recombinant Human ADP-ribosylation factor-like protein 8A protein (20-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-like protein 8A (ADP-ribosylation factor-like protein 10B) (Novel small G protein indispensable for equal chromosome segregation 2), ARL8A, ARL10B GIE2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ARL8A