



SPRR2A Antibody, FITC conjugated

Product Code	CSB-PA842658LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35326
Immunogen	Recombinant Human Small proline-rich protein 2A protein (1-72AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small proline-rich protein 2A (SPR-2A) (2-1), SPRR2A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SPRR2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.