



SPDL1 Antibody, Biotin conjugated

Product Code	CSB-PA842662LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EA4
Immunogen	Recombinant Human Protein Spindly protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Spindly (hSpindly) (Arsenite-related gene 1 protein) (Coiled-coil domain-containing protein 99) (Rhabdomyosarcoma antigen MU-RMS-40.4A) (Spindle apparatus coiled-coil domain-containing protein 1), SPDL1, CCDC99
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPDL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.