





WDR34 Antibody

Product Code	CSB-PA842668LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EX3
Immunogen	Recombinant Human WD repeat-containing protein 34 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 34, WDR34
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WDR34

Image



Western blot
 All lanes: WDR34 antibody at 1.5µg/ml
 Lane 1: PC-3 whole cell lysate
 Lane 2: MDA-MB-231 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 58 kDa
 Observed band size: 58 kDa



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA842668LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.