



WDR34 Antibody, HRP conjugated

Product Code	CSB-PA842668LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EX3
Immunogen	Recombinant Human WD repeat-containing protein 34 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 34, WDR34
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WDR34
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.