



SPG11 Antibody, Biotin conjugated

Product Code	CSB-PA842692LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JI7
Immunogen	Recombinant Human Spatacsin protein (2019-2318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spatacsin (Colorectal carcinoma-associated protein) (Spastic paraplegia 11 protein), SPG11, KIAA1840
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPG11