



ZDHHC15 Antibody, FITC conjugated

Product Code	CSB-PA842720LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MV8
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC15 protein (239-337AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC15 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 15) (DHHC-15), ZDHHC15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZDHHC15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.