



PRELID3A Antibody

Product Code	CSB-PA842722LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96N28
Immunogen	Recombinant Human PRELI domain containing protein 3A protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRELI domain containing protein 3A (Protein slowmo homolog 1), PRELID3A, C18orf43 SLMO1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRELID3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.