





RASL10B Antibody

Product Code	CSB-PA842758LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96S79
Immunogen	Recombinant Human Ras-like protein family member 10B protein (15-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-like protein family member 10B (Ras-like protein VTS58635) (Ras-related protein 17) (RRP17), RASL10B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASL10B

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA842758LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA842758LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.