





PPP1R16B Antibody

Product Code	CSB-PA842762LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96T49
Immunogen	Recombinant Human Protein phosphatase 1 regulatory inhibitor subunit 16B protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1 regulatory inhibitor subunit 16B (Ankyrin repeat domain-containing protein 4) (CAAX box protein TIMAP) (TGF-beta-inhibited membrane-associated protein) (hTIMAP), PPP1R16B, ANKRD4 KIAA0823
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP1R16B

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA842762LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA842762LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.