



ZFP82 Antibody, HRP conjugated

Product Code	CSB-PA843135LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N141
Immunogen	Recombinant Human Zinc finger protein 82 homolog protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 82 homolog (Zfp-82) (Zinc finger protein 545), ZFP82, KIAA1948 ZNF545
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZFP82
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.