



MINDY1 Antibody, HRP conjugated

Product Code	CSB-PA843272LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5J2
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase MINDY-1 protein (3-278AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase MINDY-1 (EC 3.4.19.12) (Deubiquitinating enzyme MINDY-1) (Protein FAM63A), MINDY1, FAM63A KIAA1390
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MINDY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.