



PAF1 Antibody, HRP conjugated

Product Code	CSB-PA843285LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N7H5
Immunogen	Recombinant Human RNA polymerase II-associated factor 1 homolog protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA polymerase II-associated factor 1 homolog (hPAF1) (Pancreatic differentiation protein 2), PAF1, PD2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PAF1