



# PNPLA1 Antibody, Biotin conjugated

Product Code	CSB-PA843295LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8W4
Immunogen	Recombinant Human Patatin-like phospholipase domain-containing protein 1 protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Patatin-like phospholipase domain-containing protein 1 (EC 3.1.1.-), PNPLA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	PNPLA1