




LRRC57 Antibody

Product Code	CSB-PA843302LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9N7
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 57 protein (1-239AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 57, LRRC57
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRRC57

Image



Western blot
 All lanes: LRRC57 antibody at 4μg/ml
 Lane 1: Mouse liver tissue
 Lane 2: Mouse kidney tissue
 Lane 3: Mouse lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 27 kDa
 Observed band size: 27 kDa