



LRRC57 Antibody, HRP conjugated

Product Code	CSB-PA843302LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9N7
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 57 protein (1-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 57, LRRC57
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRRC57
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.