



C1orf64 Antibody

Product Code	CSB-PA843340LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEQ6
Immunogen	Recombinant Human Uncharacterized protein C1orf64 protein (1-169AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Steroid receptor-associated and regulated protein (ER-related factor) (Steroid receptor-regulated protein), SRARP, C1orf64 ERRF SSPR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SRARP