



C1orf64 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA843340LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8NEQ6 |
| Immunogen | Recombinant Human Uncharacterized protein C1orf64 protein (1-169AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Steroid receptor-associated and regulated protein (ER-related factor) (Steroid receptor-regulated protein), SRARP, C1orf64 ERRF SSPR |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | SRARP |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |