




NUP43 Antibody

Product Code	CSB-PA844045LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFH3
Immunogen	Recombinant Human Nucleoporin Nup43 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoporin Nup43 (Nup107-160 subcomplex subunit Nup43) (p42), NUP43
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NUP43

Image



Western blot

All lanes: NUP43 antibody at 4µg/ml

Lane 1: A431 whole cell lysate

Lane 2: HL60 whole cell lysate

Lane 3: THP-1 whole cell lysate


Lane 4: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 43, 32 kDa

Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded

human testis tissue using CSB-

PA844045LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.