



GPR156 Antibody, Biotin conjugated

Product Code	CSB-PA844046LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFN8
Immunogen	Recombinant Human Probable G-protein coupled receptor 156 protein (614-814AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 156 (G-protein coupled receptor PGR28) (GABAB-related G-protein coupled receptor), GPR156, GABABL PGR28
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR156
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.