



# GPR156 Antibody, Biotin conjugated

Product Code	CSB-PA844046LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFN8
Immunogen	Recombinant Human Probable G-protein coupled receptor 156 protein (614-814AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 156 (G-protein coupled receptor PGR28) (GABAB-related G-protein coupled receptor), GPR156, GABABL PGR28
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR156