



MFSD8 Antibody, HRP conjugated

Product Code	CSB-PA844087LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHS3
Immunogen	Recombinant Human Major facilitator superfamily domain-containing protein 8 protein (1-40AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major facilitator superfamily domain-containing protein 8 (Ceroid-lipofuscinosis neuronal protein 7), MFSD8, CLN7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MFSD8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.