



MFSD8 Antibody, HRP conjugated

Product Code	CSB-PA844087LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHS3
Immunogen	Recombinant Human Major facilitator superfamily domain-containing protein 8 protein (1-40AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major facilitator superfamily domain-containing protein 8 (Ceroid-lipofuscinosis neuronal protein 7), MFSD8, CLN7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MFSD8