




NFKBID Antibody

Product Code	CSB-PA844092LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NI38
Immunogen	Recombinant Human NF-kappa-B inhibitor delta protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NF-kappa-B inhibitor delta (I-kappa-B-delta) (Ikb-delta) (IkappaBdelta) (IkappaBNS) (T-cell activation NFKB-like protein) (TA-NFKBH), NFKBID, IKBNS
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NFKBID

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, MCF-7 whole cell lysate

All lanes: NFKBID antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 34, 51, 16 kDa

Observed band size: 34 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.