



UBTD2 Antibody

Product Code	CSB-PA845144LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUN7
Immunogen	Recombinant Human Ubiquitin domain-containing protein 2 protein (45-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin domain-containing protein 2 (Dendritic cell-derived ubiquitin-like protein) (DC-UbP) (Ubiquitin-like protein SB72), UBTD2, DCUBP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UBTD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.