



Scgb3a2 Antibody, Biotin conjugated

Product Code	CSB-PA846028LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q920H1
Immunogen	Recombinant Mouse Secretoglobin family 3A member 2 protein (22-139AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretoglobin family 3A member 2 (Pneumo secretory protein 1) (PnSP-1) (Uteroglobin-related protein 1), Scgb3a2, Pnsp1 Ugrp1
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Scgb3a2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.