



RIT1 Antibody, Biotin conjugated

Product Code	CSB-PA846116LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92963
Immunogen	Recombinant Human GTP-binding protein Rit1 protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTP-binding protein Rit1 (Ras-like protein expressed in many tissues) (Ras-like without CAAX protein 1), RIT1, RIBB RIT ROC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RIT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.