



PTP4A1 Antibody, HRP conjugated

Product Code	CSB-PA846294LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93096
Immunogen	Recombinant Human Protein tyrosine phosphatase type IVA 1 protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein tyrosine phosphatase type IVA 1 (EC 3.1.3.48) (PTP(CAAXI)) (Protein-tyrosine phosphatase 4a1) (Protein-tyrosine phosphatase of regenerating liver 1) (PRL-1), PTP4A1, PRL1 PTPCAAX1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTP4A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.