



SCMH1 Antibody, Biotin conjugated

Product Code	CSB-PA846625LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GD3
Immunogen	Recombinant Human Polycomb protein SCMH1 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycomb protein SCMH1 (Sex comb on midleg homolog 1), SCMH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SCMH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.