



PIGS Antibody, FITC conjugated

Product Code	CSB-PA847222LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96S52
Immunogen	Recombinant Human GPI transamidase component PIG-S protein (40-300AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPI transamidase component PIG-S (Phosphatidylinositol-glycan biosynthesis class S protein), PIGS
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PIGS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.