



# INTS14 Antibody, HRP conjugated

Product Code	CSB-PA847227LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96SY0
Immunogen	Recombinant Human Integrator complex subunit 14 protein (231-518AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrator complex subunit 14 (von Willebrand factor A domain-containing protein 9), INTS14, C15orf44 VWA9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	INTS14