



INTS14 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA847227LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96SY0 |
| Immunogen | Recombinant Human Integrator complex subunit 14 protein (231-518AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Integrator complex subunit 14 (von Willebrand factor A domain-containing protein 9), INTS14, C15orf44 VWA9 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | INTS14 |