



INTS14 Antibody, Biotin conjugated

Product Code	CSB-PA847227LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96SY0
Immunogen	Recombinant Human Integrator complex subunit 14 protein (231-518AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrator complex subunit 14 (von Willebrand factor A domain-containing protein 9), INTS14, C15orf44 VWA9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	INTS14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.