



C1orf87 Antibody

Product Code	CSB-PA847603LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0U7
Immunogen	Recombinant Human Uncharacterized protein C1orf87 protein (1-138AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C1orf87, C1orf87
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C1orf87
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.