



SYNE4 Antibody, Biotin conjugated

Product Code	CSB-PA847615LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N205
Immunogen	Recombinant Human Nesprin-4 protein (1-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nesprin-4 (KASH domain-containing protein 4) (KASH4) (Nuclear envelope spectrin repeat protein 4), SYNE4, C19orf46
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SYNE4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.